




# General Exam Scope

## The Only Video Scope You'll Ever Need


We are proud to introduce a Medical Video Scope system that is so versatile we are calling it the General Exam Scope. It offers superior quality and advanced functionality not found in other video systems. It comes with a Camera and Illumination System and a trigger hand piece that is so versatile it can be used in a variety of medical specialties.

This camera system has **S-Video, Composite Video, XGA outputs** and has the ability to capture four images inside of the camera itself. These captured images can be looked at one at a time or tiled in a quad fashion. This is extremely helpful during a video consultation. It is also a single cable system; both the fiber optic light and video signal are contained in the same cable.

With accessories, the General Exam Scope can be used as an Otoscope, Dermoscope, Colposcope, and much more.



• 1 Display



• 4 Display

### Basic



• Side illumination

• Tongue depressor & colposc



### ACCESSORIES & SOFTWARE

#### DCS-104T

X20 IRIS CAP



SKIN



EAR DRUM, TURBINATE



X35 IRIS CAP



#### Option

LARYNX (OPTION)



ANUS, BOWELS (OPTION)



TONSIL, CERVIX



Distributed Exclusively by Second Opinion Telemedicine Solutions, Inc.

310-802-6300 email: sosales@2opinion.com www.2opinion.com

3848 Del Amo Blvd, Ste 303, Torrance, CA 90503